

Product Datasheet

Rat C10orf59 ELISA Kit (orb411141)

Description	Rat Renalase/ Monoamine Oxidase-C(RNLS/ MAO-C)ELISA Kit
Reactivity	Rat
Range	62.5 ng/mL-4000 ng/mL
Target	C10orf59
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q5U2W9
Sensitivity	15.6 ng/mL
Expiration Date	Please enquire.