

Product Datasheet

Rat CXCL1 ELISA Kit (orb411150)

Description	Rat growth-regulated oncogene α /melanoma growth stimulating activity,GRO α /MGSA Elisa Kit
Reactivity	Rat
Range	125 pg/mL-8000 pg/mL
Target	CXCL1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P14095
Sensitivity	31.25 pg/mL
Expiration Date	Please enquire.