

Product Datasheet

Rat Glucocorticoid Receptor ELISA Kit (orb411185)

Description	Rat glucocorticoid receptor,GR ELISA Kit
Reactivity	Rat
Range	15.6 pg/mL-1000 pg/mL
Target	NR3C1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P06536
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.