

## Product Datasheet

### Rat Dehydroepiandrosterone ELISA Kit (orb411201)

|                          |                                                 |
|--------------------------|-------------------------------------------------|
| <b>Description</b>       | Rat<br>Dehydroepiandrosterone,DHEA<br>ELISA Kit |
| <b>Reactivity</b>        | Rat                                             |
| <b>Range</b>             | 0.05 ng/mL-10 ng/mL                             |
| <b>Target</b>            | DHEA                                            |
| <b>Note</b>              | For research use only                           |
| <b>Application notes</b> | Detection Wavelength: 450 nm                    |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates            |
| <b>Sensitivity</b>       | 0.025 ng/mL                                     |
| <b>Expiration Date</b>   | Please enquire.                                 |