

## Product Datasheet

### Rat RARRES2 ELISA Kit (orb411209)

|                          |                                                                   |
|--------------------------|-------------------------------------------------------------------|
| <b>Description</b>       | Rat Retinoic acid receptor responder protein 2(RARRES2) ELISA kit |
| <b>Reactivity</b>        | Rat                                                               |
| <b>Range</b>             | 0.9 ng/mL-60 ng/mL                                                |
| <b>Target</b>            | RARRES2                                                           |
| <b>Note</b>              | For research use only                                             |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                      |
| <b>Sample Types</b>      | serum, plasma, cell culture supernates, tissue homogenates        |
| <b>Sensitivity</b>       | 0.225 ng/mL                                                       |
| <b>Expiration Date</b>   | Please enquire.                                                   |