

Product Datasheet

Rat CK 18-M65 ELISA Kit (orb411221)

Description	Rat cytokeratin 18-M65(CK 18-M65) ELISA kit
Reactivity	Rat
Range	0.94 ng/mL-60 ng/mL
Target	KRT18
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma,tissue homogenates
Uniprot ID	Q5BJY9
Sensitivity	0.23 ng/mL
Expiration Date	Please enquire.