

Product Datasheet

Rat CYP7A1 ELISA Kit (orb411225)

Description	Rat cytochrome P450, family 7, subfamily A, polypeptide 1 (CYP7A1) ELISA kit
Reactivity	Rat
Range	93.75 pg/mL-6000 pg/mL
Target	CYP7A1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P18125
Sensitivity	23.4 pg/mL
Expiration Date	Please enquire.