

Product Datasheet

Rat ACE2 ELISA Kit (orb411261)

Description	Rat Angiotensin converting enzyme 2, ACE2 ELISA Kit
Reactivity	Rat
Range	0.02 ng/ml - 20 ng/ml
Target	ACE2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q5EGZ1
Sensitivity	0.02 ng/ml
Expiration Date	Please enquire.