

Product Datasheet

Rat FLT1 ELISA Kit (orb411266)

Description	Rat Vascular endothelial cell growth factor receptor 1, VEGFR-1/Flt1 ELISA Kit
Reactivity	Rat
Range	1.56 pg/mL-100 pg/mL
Target	VEGFR-1/Flt1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P53767
Sensitivity	0.39 pg/mL
Expiration Date	Please enquire.