

Product Datasheet

Rat KDR ELISA Kit (orb411267)

Description	Rat Vascular endothelial cell growth factor receptor 2, VEGFR-2/Flk-1 ELISA kit
Reactivity	Rat
Range	1.56 ng/mL-100 ng/mL
Target	VEGFR-2/Flk-1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	008775
Sensitivity	0.39 ng/mL
Expiration Date	Please enquire.