

Product Datasheet

Rat HTR2A ELISA Kit (orb411282)

Description	Rat Serotonin Receptor 2A/5-hydroxytryptamine Receptor 2A(5-HT-2A)ELISA Kit
Reactivity	Rat
Range	31.25 pg/mL-2000 pg/mL
Target	HTR2A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P14842
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.