

Product Datasheet

Rat GPIIbIIIa/CD41+CD61 ELISA Kit (orb411289)

Description	Rat platelet membrane glycoprotein IIbIIIa,GP-IIbIIIa ELISA Kit
Reactivity	Rat
Range	6.25 ng/mL-400 ng/mL
Target	GPIIbIIIa/CD41+CD61
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Sensitivity	1.56 ng/mL
Expiration Date	Please enquire.