

## Product Datasheet

### Rat Circulating immune complex ELISA Kit (orb411297)

|                          |                                              |
|--------------------------|----------------------------------------------|
| <b>Description</b>       | Rat circulating immune complex(CIC)ELISA Kit |
| <b>Reactivity</b>        | Rat                                          |
| <b>Range</b>             | Request Information                          |
| <b>Target</b>            | CIC                                          |
| <b>Note</b>              | For research use only                        |
| <b>Application notes</b> | Detection Wavelength: 450 nm                 |
| <b>Sensitivity</b>       | Request Information                          |
| <b>Expiration Date</b>   | Please enquire.                              |