

Product Datasheet

Rat CHRM1 ELISA Kit (orb411298)

Description	Rat Muscarinic acetylcholine receptor M1(CHRM1) ELISA kit
Reactivity	Rat
Range	0.27 ng/ml- 200ng/ml
Target	CHRM1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates.
Uniprot ID	P08482
Sensitivity	0.27 ng/ml
Expiration Date	Please enquire.