

Product Datasheet

Rat IGF1R ELISA Kit (orb411303)

Description	Rat Insulin-like Growth Factors 1 Receptor,IGF-1R ELISA Kit
Reactivity	Rat
Range	3.12 pg/mL-200 pg/mL
Target	IGF1R
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P24062
Sensitivity	0.8 pg/mL
Expiration Date	Please enquire.