

Product Datasheet

Rat IGFBP1 ELISA Kit (orb411306)

Description	Rat Insulin-like growth factor binding protein 1,IGFBP-1 ELISA Kit
Reactivity	Rat
Range	15.6 pg/mL-1000 pg/mL
Target	IGFBP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P21743
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.