

Product Datasheet

Rat IGFBP2 ELISA Kit (orb411307)

Description	Rat Insulin-like growth factor binding protein 2,IGFBP-2 ELISA kit
Reactivity	Rat
Range	0.625 ng/mL-40 ng/mL
Target	IGFBP2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P12843
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.