

Product Datasheet

Rat IGFBP3 ELISA Kit (orb411308)

Description	Rat Insulin-like growth factor binding protein 3,IGFBP-3 ELISA kit
Reactivity	Rat
Range	2.7 ng/ml - 2000 ng/ml
Target	IGFBP3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P15473
Sensitivity	2.7 ng/ml
Expiration Date	Please enquire.