

Product Datasheet

Rat IDO ELISA Kit (orb411336)

Description	Rat Indoleamine 2,3-dioxygenase,IDO ELISA Kit
Reactivity	Rat
Range	31.25 pg/mL-2000 pg/mL
Target	IDO1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9ERD9
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.