

Product Datasheet

Rat Lp-PLA2 ELISA Kit (orb411357)

Description	Rat lipoprotein-associated phospholipase A2,Lp-PLA2 ELISA Kit
Reactivity	Rat
Range	0.78 IU/mL-50 IU/mL
Target	Lp-PLA2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma
Uniprot ID	Q5M7T7
Sensitivity	0.195 IU/mL
Expiration Date	Please enquire.