

Product Datasheet

Rat KISS1 ELISA Kit (orb411375)

Description	Rat Kisspeptin 1(KISS1)ELISA Kit
Reactivity	Rat
Range	0.156 ng/mL-10 ng/mL
Target	KISS1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q7TSB7
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.