

Product Datasheet

Rat TGF beta 1 ELISA Kit (orb411378)

Description	Rat Transforming Growth Factor β 1(TGF- β 1) ELISA kit
Reactivity	Rat
Range	6.25 pg/mL-400 pg/mL
Target	TGFB1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P17246
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.