

## Product Datasheet

### Sheep ACTH ELISA Kit (orb411455)

|                          |                                                     |
|--------------------------|-----------------------------------------------------|
| <b>Description</b>       | Sheep adrenocorticotrophic hormone (ACTH) ELISA kit |
| <b>Reactivity</b>        | Sheep                                               |
| <b>Range</b>             | 0.097 pg/mL-400 pg/mL                               |
| <b>Target</b>            | ACTH                                                |
| <b>Note</b>              | For research use only                               |
| <b>Application notes</b> | Detection Wavelength: 450 nm                        |
| <b>Sample Types</b>      | serum, plasma                                       |
| <b>Sensitivity</b>       | 0.097 pg/mL                                         |
| <b>Expiration Date</b>   | Please enquire.                                     |