

Product Datasheet

Sheep Prolactin ELISA Kit (orb411459)

Description	Sheep Prolactin/Luteotropic Hormone(PRL/LTH) ELISA Kit
Reactivity	Sheep
Range	5 µIU/mL-2000 µIU/mL
Target	PRL
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P01240
Sensitivity	5 µIU/mL
Expiration Date	Please enquire.