

## Product Datasheet

### Sheep GHRL ELISA Kit (orb411460)

|                          |                                                        |
|--------------------------|--------------------------------------------------------|
| <b>Description</b>       | Sheep ghrelin/obestatin prepropeptide (GHRL) ELISA kit |
| <b>Reactivity</b>        | Sheep                                                  |
| <b>Range</b>             | 62.5 pg/mL-4000 pg/mL                                  |
| <b>Target</b>            | GHRL                                                   |
| <b>Note</b>              | For research use only                                  |
| <b>Application notes</b> | Detection Wavelength: 450 nm                           |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                      |
| <b>Sensitivity</b>       | 15.6 pg/mL                                             |
| <b>Expiration Date</b>   | Please enquire.                                        |