

## Product Datasheet

### Sheep TNF ELISA Kit (orb411490)

**Description** Sheep Tumor Necrosis Factor  $\alpha$ , TNF- $\alpha$   
ELISA Kit

**Reactivity** Sheep

**Range** 3.125 pg/mL-50 pg/mL

**Target** TNF

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, tissue homogenates

**Uniprot ID** **P23383**

**Sensitivity** 3.125 pg/mL

**Expiration Date** Please enquire.