

Product Datasheet

HLA class II antibody (orb411560)

Description	Mouse monoclonal antibody to HLA class II
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Target	HLA class II
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody has been show to work by flow-cytometry and western blot. Dechant et al (PMID: 12393717) used a human IgA2-chimeric version of this antibody to study IgA-mediated cell killing.
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.