

Product Datasheet

Alpha-Tubulin antibody (orb411569)

Description	Human monoclonal antibody to Alpha-Tubulin
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Human MBP (microtubule-binding protein).
Target	Alpha-Tubulin
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody binds alpha-tubulin. Tubulin is the major constituent of microtubules. The alpha-chain has a non-exchangeable GTP-binding site.
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q71U36
Expiration Date	12 months from date of receipt.