

Product Datasheet

CD8 beta antibody (orb411576)

Description	Rat monoclonal antibody to CD8 beta
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC
Immunogen	Mouse CD8 (Iyt-3).
Target	CD8 beta
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody recognizes murine CD8 beta.
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P10300
Expiration Date	12 months from date of receipt.