

Product Datasheet

VEGF antibody (orb411577)

Description	Mouse monoclonal antibody to VEGF
Species/Host	Mouse
Reactivity	Human, Mouse, Primate
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Human VEGF189 expressed in E. coli.
Target	VEGF
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	VG76e will detect both native and denatured VEGF, an angiogenic growth factor and a prognostic indicator for cancer when detected in serum.
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P15692
Expiration Date	12 months from date of receipt.