

Product Datasheet

amyloid beta antibody (orb411589)

Description	Mouse monoclonal antibody to amyloid beta
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, EM, IHC, IP, WB
Immunogen	Amyloid-beta peptide.
Target	amyloid beta
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody recognises human amyloid beta and can be used to detect amyloid in ELISAs, WB and tissue sections.
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P05067
Expiration Date	12 months from date of receipt.