



Product Datasheet

GSTM2 antibody (orb411597)

Description	Rabbit polyclonal antibody to GSTM2.
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human GSTM2(NP_000839.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	27kDa
Uniprot ID	P28161
Dilution Range	WB 1:500-:2000, IHC-P: 1:50-:200, IF/ICC: 1:50:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549 cell...



Immunohistochemical staining of human co...