



Product Datasheet

PPP3R1 antibody (orb411631)

Description	Rabbit polyclonal antibody to PPP3R1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant full length protein of human PPP3R1
Target	PPP3R1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P63100, Q63810, P63098
Entrez	5534, 29748, 19058
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IP: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of PPP3R1 expressi...



Immunofluorescent analysis of PPP3R1 sta...