



Product Datasheet

RASSF6 antibody (orb41235)


Description	RASSF6 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Ras association domain-containing protein 6 protein (1-303AA)
Target	RASSF6
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6ZTQ3
Dilution Range	IHC-P:1:20-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human prostate tissue.



Immunohistochemical staining of human kidney tissue.



Immunohistochemical staining of human ovary tissue.