



Product Datasheet

KMT2E antibody (orb412434)


Description	Rabbit polyclonal antibody to KMT2E
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, IP, WB
Immunogen	Recombinant full length protein of human MLL5
Target	KMT2E
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q3UG20 , Q8IZD2
Entrez	55904 , 69188
Dilution Range	WB: 1:500-2000, IF/ICC: 1:50-200, IP: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of MLL5 expression...



Immunohistochemical analysis of MLL5 staining...



Immunofluorescent analysis of MLL5 staining...