



Product Datasheet

Histone H3 antibody (orb412573)


Description	Rabbit polyclonal antibody to Histone H3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, IP, WB
Immunogen	KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K27 of human Histone H3 protein. The exact sequence is proprietary.
Target	HIST3H3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q16695
Entrez	8290
Dilution Range	WB: 1:500-2000, IP: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone H3 (Acetylated)



Immunohistochemical analysis of Histone H3 (Acetylated)



Immunofluorescent analysis of Histone H3 (Acetylated)