


## Product Datasheet

### EDIL3 antibody (orb41296)

|                            |                                                                                                       |
|----------------------------|-------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | EDIL3 antibody                                                                                        |
| <b>Species/Host</b>        | Rabbit                                                                                                |
| <b>Reactivity</b>          | Human                                                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                                          |
| <b>Tested Applications</b> | ELISA, IHC                                                                                            |
| <b>Immunogen</b>           | Recombinant Human EGF-like repeat and discoidin I-like domain-containing protein 3 protein (60-300AA) |
| <b>Target</b>              | EDIL3                                                                                                 |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.                                                     |
| <b>Form/Appearance</b>     | Liquid                                                                                                |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                         |
| <b>Note</b>                | For research use only                                                                                 |
| <b>Isotype</b>             | IgG                                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                            |
| <b>Uniprot ID</b>          | <b>O43854</b>                                                                                         |
| <b>Dilution Range</b>      | IHC-P:1:20-200                                                                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                       |



Immunohistochemistry  
of paraffin-embedde...