



Product Datasheet

COX IV/COX4I1 Antibody (orb412995)


Description	COX IV/COX4I1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, IHC, WB
Immunogen	E. coli-derived human COX IV recombinant protein (Position: Q59-K169).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	17 kDa
Uniprot ID	P13073
Dilution Range	Western blot, 0.1-0.5 µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1 µg/ml Immunocytochemistry/Immunofluorescence, 2 µg/ml Flow Cytometry, 1-3 µg/1x10 ⁶ cells Direct ELISA, 0.1-0.5 µg/ml
Expiration Date	12 months from date of receipt.



Flow Cytometry analysis of U937 cells US...



WB analysis of COX IV using anti-COX IV ...



WB analysis of COX IV using anti-COX IV ...