

Product Datasheet

SNRPB2 antibody (orb41302)

Description	SNRPB2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human U2 small nuclear ribonucleoprotein B\\'\\'\\' protein (1-150AA)
Target	SNRPB2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08579
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.

