

Product Datasheet

Cannabinoid Receptor 1/CNR1 Antibody (orb413041)

Description

Cannabinoid Receptor 1/CNR1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, ICC, IHC, WB

Applications

Immunogen

E. coli-derived human Cannabinoid Receptor 1 recombinant protein (Position: M1-Q75).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

60 kDa

Uniprot ID


P21554

Dilution Range

Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml Immunohistochemistry (Frozen Section), 0.5-1µg/ml Immunocytochemistry, 0.5-1µg/ml Flow Cytometry, 1-3µg/1x10⁶ cells Direct ELISA, 0.1-0.5µg/ml

Expiration Date

12 months from date of receipt.




Flow Cytometry analysis of HeLa cells

US...



WB analysis of Cannabinoid Receptor 1

US...



IHC analysis of Cannabinoid Receptor 1

U...