

Product Datasheet

Insulin degrading enzyme/IDE Antibody (orb413045)

Description

Insulin degrading enzyme/IDE Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, ICC, IHC, WB

Applications

Immunogen

E. coli-derived human IDE recombinant protein (Position: F485-K756).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

118 kDa

Uniprot ID

P14735


Dilution Range

Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml Immunohistochemistry (Frozen Section), 0.5-1µg/ml Immunocytochemistry, 0.5-1µg/ml Flow Cytometry, 1-3µg/1x10⁶ cells Direct ELISA, 0.1-0.5µg/ml


Expiration Date

12 months from date of receipt.


1 2 3 4' 5 6 7 8 9 10



Western blot analysis of human Hela whole...



human placenta tissue lysates(Lane 2), h...



IHC analysis of human mammary cancer tis...