



Product Datasheet

NFIB/NF1B2 Antibody (orb413065)


Description	NFIB/NF1B2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence of human NFIB/NF1B2 (ELVRVSRTPITQGTGVNFPIGEIPSQPYHDMNSGVNLQR).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	50 kDa
Uniprot ID	O00712
Dilution Range	Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml
Expiration Date	12 months from date of receipt.



WB analysis of NFIB/NF1B2 using anti-NFIB/NF1B2 antibody.



IHC analysis of NFIB/NF1B2 using anti-NFIB/NF1B2 antibody.



IHC analysis of NFIB/NF1B2 using anti-NFIB/NF1B2 antibody.