


Product Datasheet

Human PLA2G7/Lp-PLA2 ELISA Kit (orb413105)

Description	Human PLA2G7/Lp-PLA2 ELISA Kit
Reactivity	Human
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: F22-N441
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q13093
Sensitivity	20 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is 16-200 ng/ml.
Expiration Date	Please enquire.



Human
PLA2G7
PicoKine
ELISA Kit
standard...