

Product Datasheet

Human Hyaluronidase1/HYAL1 ELISA Kit (orb413106)

Description Human Hyaluronidase1/HYAL1 ELISA Kit

Reactivity Human

Range 62.5 pg/ml - 4,000 pg/ml

Tested Applications ELISA

Immunogen Expression system for standard: NS0; Immunogen sequence: F22-W435

Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

Note For research use only


Sample Types cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).

Uniprot ID [Q12794](#)

Sensitivity 10 pg/ml

Dilution Range Dilution ratio of 1:100, concentration in serum and plasma is around 100 ng/ml.

Expiration Date Please enquire.



Human
HYAL1
PicoKine
ELISA Kit
standard ...