

Product Datasheet

PJ34(free base) (orb413175)

Description	PJ34 is a potent specific inhibitor of PARP1/2 with IC50 of 110 nM and 86 nM, respectively.
Note	For research use only
Formula	$C_{17}H_{17}N_3O_2$
MW	295.34
CAS Number	[344458-19-1]
SMILES	<chem>CN(CC(NC1C=CC2=C(C3=CC=CC=C3C(N2)=O)C=1)=O)C</chem>
Expiration Date	12 months from date of receipt.