



Product Datasheet

BST1 antibody (orb41355)

Description	BST1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 2 protein (29-293AA)
Target	BST1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q10588
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis
of Mouse lung tiss...



Immunohistochemical
staining of human
te...