

Product Datasheet

Ubiquitin antibody (orb413632)

Description	Rabbit Recombinant monoclonal antibody against Ubiquitin conjugated to Biotin
Species/Host	Rabbit
Reactivity	All
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Ubiquitin conjugated lysozyme
Target	Ubiquitin
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% Sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SMC-548 was sufficient for detection of Ubiquitin in 2 µg of purified poly-ubiquitin chains (Poly-Ubiquitin: Ub1-Ub7 WT Chains K63-linked, k48-linked) by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Isotype	IgG
Clonality	Recombinant
Uniprot ID	P0CG47, P0CG48
Entrez	7314
Dilution Range	WB (1:1000); ELISA (1:1000)
Expiration Date	12 months from date of receipt.