


Product Datasheet

TNFSF15 antibody (orb41367)

Description	TNFSF15 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Tumor necrosis factor ligand superfamily member 15 protein (72-251AA)
Target	TNFSF15
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95150
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse lung tiss...