

Product Datasheet

**DYKDDDDK Tag antibody (Equivalent to Sigma's
Anti-FLAG M2 antibody) [Out of stock]
(orb413733)**

Description

Rabbit polyclonal antibody

Species/Host

Rabbit

Reactivity

All

Conjugation

ATTO-594

Tested

WB

Applications
Immunogen

Synthetic peptide corresponding to the DYKDDDDK Tag, conjugated to Keyhole Limpet Haemocyanin (KLH).

Target

DYKDDDDK Tag

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

A 1:1000 dilution of SPC-720 was sufficient for detection of DYKDDDDK Tag in 15 µg of DYKDDDDK (FLAG) Positive Control Lysate cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

Dilution Range

WB (1:1000)

Expiration Date

12 months from date of receipt.