

Product Datasheet

TRHDE antibody (orb41377)

Description	TRHDE antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Thyrotropin-releasing hormone-degrading ectoenzyme protein (62-200AA)
Target	TRHDE
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UKU6
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.